

Roth Prescription

						10個/パック		10個/パック	
上顎						V-Slot無し		V-Slot付	
部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
中切歯	+12	+5	0	3.5mm	R	36-U1-8R	36-U1-2R	38-U1-8R	38-U1-2R
					L	36-U1-8L	36-U1-2L	38-U1-8L	38-U1-2L
中切歯 ハイトルク	+22	+5	0	3.5mm	R	39-U1-8R	39-U1-2R	---	---
					L	39-U1-8L	39-U1-2L	---	---
側切歯	+8	+9	0	3.0mm	R	36-U2-8R	36-U2-2R	38-U2-8R	38-U2-2R
					L	36-U2-8L	36-U2-2L	38-U2-8L	38-U2-2L
側切歯 ハイトルク	+14	+9	0	3.0mm	R	39-U2-8R	39-U2-2R	---	---
					L	39-U2-8L	39-U2-2L	---	---
犬歯	0	+11	0	3.3mm	R	36-U3-8R	36-U3-2R	38-U3-8R	38-U3-2R
					L	36-U3-8L	36-U3-2L	38-U3-8L	38-U3-2L
犬歯 フック付					R	36-U3-8RK	36-U3-2RK	38-U3-8RK	38-U3-2RK
					L	36-U3-8LK	36-U3-2LK	38-U3-8LK	38-U3-2LK
小臼歯	-7	0	0	3.0mm	Uni	36-U4-8U	36-U4-2U	38-U4-8U	38-U4-2U
小臼歯 フック付					R	36-U4-8RK	36-U4-2RK	38-U4-8RK	38-U4-2RK
					L	36-U4-8LK	36-U4-2LK	38-U4-8LK	38-U4-2LK

下顎

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
前歯	0	0	0	2.5mm	Uni	36-L1-8U	36-L1-2U	38-L1-8U	38-L1-2U
犬歯	-11	+5	0	3.0mm	R	36-L3-8R	36-L3-2R	38-L3-8R	38-L3-2R
					L	36-L3-8L	36-L3-2L	38-L3-8L	38-L3-2L
犬歯 フック付					R	36-L3-8RK	36-L3-2RK	38-L3-8RK	38-L3-2RK
					L	36-L3-8LK	36-L3-2LK	38-L3-8LK	38-L3-2LK
第1小臼歯	-17	0	0	3.0mm	Uni	36-L4-8U	36-L4-2U	38-L4-8U	38-L4-2U
第1小臼歯 フック付					R	36-L4-8RK	36-L4-2RK	38-L4-8RK	38-L4-2RK
					L	36-L4-8LK	36-L4-2LK	38-L4-8LK	38-L4-2LK
第2小臼歯	-22	0	0	3.0mm	Uni	36-L5-8U	36-L5-2U	38-L5-8U	38-L5-2U
第2小臼歯 フック付					R	36-L5-8RK	36-L5-2RK	38-L5-8RK	38-L5-2RK
					L	36-L5-8LK	36-L5-2LK	38-L5-8LK	38-L5-2LK

症例セット(20個/セット)

内容	V-Slot無し		V-Slot付	
	.018	.022	.018	.022
上下顎5~5 フックなし	36K-NHK-18-01	36K-NHK-22-01	38K-NHK-18	38K-NHK-22
上下顎5~5 3のみフック付	36K-3-18-01	36K-3-22-01	38K-3-18	38K-3-22
上下顎5~5 4, 5のみフック付	36K-45-18-01	36K-45-22-01	38K-45-18	38K-45-22
上下顎5~5 3, 4, 5フック付	36K-345-18-01	36K-345-227-01	38K-345-18	38K-345-22

* 全てのトルク、オフセット等の角度の単位は「°」です。

医療機器認証番号: 222AKBZX00123000

miniPrevail™ BRACKET SYSTEM



MBT Prescription

上顎

10個/パック

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022
中切歯	+17	+4	0	3.9mm	R	37-U1-8R	37-U1-2R
					L	37-U1-8L	37-U1-2L
側切歯	+10	+8	0	3.3mm	R	37-U2-8R	37-U2-2R
					L	37-U2-8L	37-U2-2L
犬歯	-7	+8	0	3.3mm	R	37-U3-8R	37-U3-2R
					L	37-U3-8L	37-U3-2L
犬歯 フック付					R	37-U3-8RK	37-U3-2RK
					L	37-U3-8LK	37-U3-2LK
犬歯	0	+8	0	3.3mm	R	37-U3-8R0	37-U3-2R0
					L	37-U3-8L0	37-U3-2L0
犬歯 フック付					R	37-U3-8RK0	37-U3-2RK0
					L	37-U3-8LK0	37-U3-2LK0
小臼歯	-7	0	0	3.5mm	Uni	37-U4-8U	37-U4-2U
小臼歯 フック付					R	37-U4-8RK	37-U4-2RK
					L	37-U4-8LK	37-U4-2LK



下顎

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022
前歯	-6	0	0	2.6mm	Uni	37-L1-8U	37-L1-2U
犬歯	-6	+3	0	3.5mm	R	37-L3-8R	37-L3-2R
					L	37-L3-8L	37-L3-2L
犬歯 フック付					R	37-L3-8RK	37-L3-2RK
					L	37-L3-8LK	37-L3-2LK
犬歯	0	+3	0	3.5mm	R	37-L3-8R0	37-L3-2R0
					L	37-L3-8L0	37-L3-2L0
犬歯 フック付					R	37-L3-8RK0	37-L3-2RK0
					L	37-L3-8LK0	37-L3-2LK0
第1小臼歯	-12	0	0	3.5mm	Uni	37-L4-8U	37-L4-2U
第1小臼歯 フック付					R	37-L4-8RK	37-L4-2RK
					L	37-L4-8LK	37-L4-2LK
第2小臼歯	-17	0	0	3.5mm	Uni	37-L5-8U	37-L5-2U
第2小臼歯 フック付					R	37-L5-8RK	37-L5-2RK
					L	37-L5-8LK	37-L5-2LK



症例セット(20個/セット)

内容	.018	.022
上下顎5~5 フックなし	18K-NHK-18-01	18K-NHK-22-01
上下顎5~5 3のみフック付	18K-3-18-01	18K-3-22-01
上下顎5~5 4, 5のみフック付	18K-45-18-01	18K-45-22-01
上下顎5~5 3, 4, 5フック付	18K-345-18-01	18K-345-22-01

ComfortZone™

LOW FRICTION

by Orthodontic Design and Production, Inc.



Roth Prescription

部位	Torque	Angle	Rotation	M/D	R/L	アンカーロックベース 10個/パック		メッシュベース 10個/パック	
						.018	.022	.018	.022
中切歯	+12	+5	0	3.5mm	R	16-U1-8R	16-U1-2R	26-U1-8R	26-U1-2R
					L	16-U1-8L	16-U1-2L	26-U1-8L	26-U1-2L
中切歯 ハイトルク	+22	+5	0	3.5mm	R	17-U1-8R	17-U1-2R	---	---
					L	17-U1-8L	17-U1-2L	---	---
側切歯	+8	+9	0	3.0mm	R	16-U2-8R	16-U2-2R	26-U2-8R	26-U2-2R
					L	16-U2-8L	16-U2-2L	26-U2-8L	26-U2-2L
側切歯 ハイトルク	+14	+9	0	3.0mm	R	17-U2-8R	17-U2-2R	---	---
					L	17-U2-8L	17-U2-2L	---	---
犬歯	0	+11	0	3.3mm	R	16-U3-8R	16-U3-2R	26-U3-8R	26-U3-2R
					L	16-U3-8L	16-U3-2L	26-U3-8L	26-U3-2L
犬歯 フック付					R	16-U3-8RK	16-U3-2RK	26-U3-8RK	26-U3-2RK
					L	16-U3-8LK	16-U3-2LK	26-U3-8LK	26-U3-2LK
小臼歯	-7	0	0	3.0mm	Uni	16-U4-8U	16-U4-2U	26-U4-8U	26-U4-2U
小臼歯 フック付					R	16-U4-8RK	16-U4-2RK	26-U4-8RK	26-U4-2RK
					L	16-U4-8LK	16-U4-2LK	26-U4-8LK	26-U4-2LK
下顎									
部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
前歯	0	0	0	2.5mm	Uni	16-L1-8U	16-L1-2U	26-L1-8U	26-L1-2U
犬歯	-11	+5	0	3.0mm	R	16-L3-8R	16-L3-2R	26-L3-8R	26-L3-2R
					L	16-L3-8L	16-L3-2L	26-L3-8L	26-L3-2L
犬歯 フック付					R	16-L3-8RK	16-L3-2RK	26-L3-8RK	26-L3-2RK
					L	16-L3-8LK	16-L3-2LK	26-L3-8LK	26-L3-2LK
第1小臼歯	-17	0	0	3.0mm	Uni	16-L4-8U	16-L4-2U	26-L4-8U	26-L4-2U
第1小臼歯 フック付					R	16-L4-8RK	16-L4-2RK	26-L4-8RK	26-L4-2RK
					L	16-L4-8LK	16-L4-2LK	26-L4-8LK	26-L4-2LK
第2小臼歯	-22	0	0	3.0mm	Uni	16-L5-8U	16-L5-2U	26-L5-8U	26-L5-2U
第2小臼歯 フック付					R	16-L5-8RK	16-L5-2RK	26-L5-8RK	26-L5-2RK
					L	16-L5-8LK	16-L5-2LK	26-L5-8LK	26-L5-2LK

症例セット(20個/セット)

内容	アンカーロックベース		メッシュベース	
	.018	.022	.018	.022
上下顎5~5 フックなし	16K-NHK-18	16K-NHK-22	26K-NHK-18	26K-NHK-22
上下顎5~5 3のみフック付	16K-3-18	16K-3-22	26K-3-18	26K-3-22
上下顎5~5 4, 5のみフック付	16K-45-18	16K-45-22	26K-45-18	26K-45-22
上下顎5~5 3, 4, 5フック付	16K-345-18	16K-345-22	26K-345-18	26K-345-22

* 全てのトルク、オフセット等の角度の単位は「°」です。

医療機器認証番号: 222AKBZX00123000



10個/パック



10個/パック

MBT Prescription

上顎						アンカーロックベース		メッシュベース	
部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
中切歯	+17	+4	0	3.9mm	R	18-U1-8R	18-U1-2R	28-U1-8R	28-U1-2R
					L	18-U1-8L	18-U1-2L	28-U1-8L	28-U1-2L
側切歯	+10	+8	0	3.3mm	R	18-U2-8R	18-U2-2R	28-U2-8R	28-U2-2R
					L	18-U2-8L	18-U2-2L	28-U2-8L	28-U2-2L
犬歯	-7	+8	0	3.3mm	R	18-U3-8R	18-U3-2R	28-U3-8R	28-U3-2R
					L	18-U3-8L	18-U3-2L	28-U3-8L	28-U3-2L
犬歯 フック付					R	18-U3-8RK	18-U3-2RK	28-U3-8RK	28-U3-2RK
					L	18-U3-8LK	18-U3-2LK	28-U3-8LK	28-U3-2LK
犬歯	0	+8	0	3.3mm	R	18-U3-8R0	18-U3-2R0	28-U3-8R0	28-U3-2R0
					L	18-U3-8L0	18-U3-2L0	28-U3-8L0	28-U3-2L0
犬歯 フック付					R	18-U3-8RK0	18-U3-2RK0	28-U3-8RK0	28-U3-2RK0
					L	18-U3-8LK0	18-U3-2LK0	28-U3-8LK0	28-U3-2LK0
小臼歯	-7	0	0	3.5mm	Uni	18-U4-8U	18-U4-2U	28-U4-8U	28-U4-2U
小臼歯 フック付					R	18-U4-8RK	18-U4-2RK	28-U4-8RK	28-U4-2RK
					L	18-U4-8LK	18-U4-2LK	28-U4-8LK	28-U4-2LK
下顎						アンカーロックベース		メッシュベース	
部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
前歯	-6	0	0	2.6mm	Uni	18-L1-8U	18-L1-2U	28-L1-8U	28-L1-2U
犬歯	-6	+3	0	3.5mm	R	18-L3-8R	18-L3-2R	28-L3-8R	28-L3-2R
					L	18-L3-8L	18-L3-2L	28-L3-8L	28-L3-2L
犬歯 フック付					R	18-L3-8RK	18-L3-2RK	28-L3-8RK	28-L3-2RK
					L	18-L3-8LK	18-L3-2LK	28-L3-8LK	28-L3-2LK
犬歯	0	+3	0	3.5mm	R	18-L3-8R0	18-L3-2R0	28-L3-8R0	28-L3-2R0
					L	18-L3-8L0	18-L3-2L0	28-L3-8L0	28-L3-2L0
犬歯 フック付					R	18-L3-8RK0	18-L3-2RK0	28-L3-8RK0	28-L3-2RK0
					L	18-L3-8LK0	18-L3-2LK0	28-L3-8LK0	28-L3-2LK0
第1小臼歯	-12	0	0	3.5mm	Uni	18-L4-8U	18-L4-2U	28-L4-8U	28-L4-2U
第1小臼歯 フック付					R	18-L4-8RK	18-L4-2RK	28-L4-8RK	28-L4-2RK
					L	18-L4-8LK	18-L4-2LK	28-L4-8LK	28-L4-2LK
第2小臼歯	-17	0	0	3.5mm	Uni	18-L5-8U	18-L5-2U	28-L5-8U	28-L5-2U
第2小臼歯 フック付					R	18-L5-8RK	18-L5-2RK	28-L5-8RK	28-L5-2RK
					L	18-L5-8LK	18-L5-2LK	28-L5-8LK	28-L5-2LK

症例セット(20個/セット)

内容	アンカーロックベース		メッシュベース	
	.018	.022	.018	.022
上下顎5~5 フックなし	18K-NHK-18	18K-NHK-22	28K-NHK-18	28K-NHK-22
上下顎5~5 3のみフック付	18K-3-18	18K-3-22	28K-3-18	28K-3-22
上下顎5~5 4, 5のみフック付	18K-45-18	18K-45-22	28K-45-18	28K-45-22
上下顎5~5 3, 4, 5フック付	18K-345-18	18K-345-22	28K-345-18	28K-345-22

Edgewise Prescription

10個/パック

上顎

部位	Torque	Angle	Rotation	M/D	R/L	アンカーロックベース	
						.018	.022
中切歯	0	0	0	3.5mm	Uni	16-CEN-18	16-CEN-22
側切歯	0	0	0	3.0mm	Uni	16-LAT-18	16-LAT-22
犬歯	0	0	0	3.3mm	Uni	16-CUSP-18	16-CUSP-22
小臼歯	0	0	0	3.0mm	Uni	16-BIC-18	16-BIC-22


下顎

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022
下顎前歯	0	0	0	2.5mm	Uni	16-ANT-18	16-ANT-22

症例セット(20個/セット)

アンカーロックベース

内容	.018	.022
上下顎5~5 フックなし	16K-EDGE-18	16K-EDGE-22




*Edgewise Prescriptionは全てフック無しです。

FocusTM

BRACKET SYSTEM

by Orthodontic Design and Production, Inc.



Roth Prescription

上顎

部位	Torque	Angle	Rotation	M/D	R/L	アンカーロックベース 10個/パック		メッシュベース 10個/パック	
						.018	.022	.018	.022
中切歯	+12	+5	0	3.9mm	R	19-U1-8R	19-U1-2R	23-U1-8R	23-U1-2R
					L	19-U1-8L	19-U1-2L	23-U1-8L	23-U1-2L
側切歯	+8	+9	0	3.3mm	R	19-U2-8R	19-U2-2R	23-U2-8R	23-U2-2R
					L	19-U2-8L	19-U2-2L	23-U2-8L	23-U2-2L
犬歯	0	+11	0	3.4mm	R	19-U3-8R	19-U3-2R	23-U3-8R	23-U3-2R
					L	19-U3-8L	19-U3-2L	23-U3-8L	23-U3-2L
犬歯 フック付					R	19-U3-8RK	19-U3-2RK	23-U3-8RK	23-U3-2RK
					L	19-U3-8LK	19-U3-2LK	23-U3-8LK	23-U3-2LK
小臼歯	-7	0	0	3.0mm	Uni	19-U4-8U	19-U4-2U	23-U4-8U	23-U4-2U
小臼歯 フック付					R	19-U4-8RK	19-U4-2RK	23-U4-8RK	23-U4-2RK
					L	19-U4-8LK	19-U4-2LK	23-U4-8LK	23-U4-2LK

下顎

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
前歯	0	0	0	2.6mm	Uni	19-L1-8U	19-L1-2U	23-L1-8U	23-L1-2U
犬歯	-11	+5	0	3.4mm	R	19-L3-8R	19-L3-2R	23-L3-8R	23-L3-2R
					L	19-L3-8L	19-L3-2L	23-L3-8L	23-L3-2L
犬歯 フック付					R	19-L3-8RK	19-L3-2RK	23-L3-8RK	23-L3-2RK
					L	19-L3-8LK	19-L3-2LK	23-L3-8LK	23-L3-2LK
第1小臼歯	-17	0	0	3.4mm	Uni	19-L4-8U	19-L4-2U	23-L4-8U	23-L4-2U
第1小臼歯 フック付					R	19-L4-8RK	19-L4-2RK	23-L4-8RK	23-L4-2RK
					L	19-L4-8LK	19-L4-2LK	23-L4-8LK	23-L4-2LK
第2小臼歯	-22	0	0	3.4mm	Uni	19-L5-8U	19-L5-2U	23-L5-8U	23-L5-2U
第2小臼歯 フック付					R	19-L5-8RK	19-L5-2RK	23-L5-8RK	23-L5-2RK
					L	19-L5-8LK	19-L5-2LK	23-L5-8LK	23-L5-2LK

症例セット(20個/セット)

内容	アンカーロックベース		メッシュベース	
	.018	.022	.018	.022
上下顎5~5 フックなし	19K-NHK-18	19K-NHK-22	23K-NHK-18	23K-NHK-22
上下顎5~5 3のみフック付	19K-3-18	19K-3-22	23K-3-18	23K-3-22
上下顎5~5 4, 5のみフック付	19K-45-18	19K-45-22	23K-45-18	23K-45-22
上下顎5~5 3, 4, 5フック付	19K-345-18	19K-345-22	23K-345-18	23K-345-22

* 全てのトルク、オフセット等の角度の単位は「°」です。

医療機器認証番号: 222AKBZX00123000

FocusTM

BRACKET SYSTEM

by Orthodontic Design and Production, Inc.



MBT Prescription

上顎						アンカーロックベース 10個/パック		メッシュベース 10個/パック	
部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
中切歯	+17	+4	0	3.9mm	R	21-U1-8R	21-U1-2R	22-U1-8R	22-U1-2R
					L	21-U1-8L	21-U1-2L	22-U1-8L	22-U1-2L
側切歯	+10	+8	0	3.3mm	R	21-U2-8R	21-U2-2R	22-U2-8R	22-U2-2R
					L	21-U2-8L	21-U2-2L	22-U2-8L	22-U2-2L
犬歯	-7	+8	0	3.4mm	R	21-U3-8R	21-U3-2R	22-U3-8R	22-U3-2R
					L	21-U3-8L	21-U3-2L	22-U3-8L	22-U3-2L
犬歯 フック付					R	21-U3-8RK	21-U3-2RK	22-U3-8RK	22-U3-2RK
					L	21-U3-8LK	21-U3-2LK	22-U3-8LK	22-U3-2LK
犬歯	0	+8	0	3.4mm	R	21-U3-8R0	21-U3-2R0	22-U3-8R0	22-U3-2R0
					L	21-U3-8L0	21-U3-2L0	22-U3-8L0	22-U3-2L0
犬歯 フック付					R	21-U3-8RK0	21-U3-2RK0	22-U3-8RK0	22-U3-2RK0
					L	21-U3-8LK0	21-U3-2LK0	22-U3-8LK0	22-U3-2LK0
小臼歯	-7	0	0	3.4mm	Uni	21-U4-8U	21-U4-2U	22-U4-8U	22-U4-2U
小臼歯 フック付					R	21-U4-8RK	21-U4-2RK	22-U4-8RK	22-U4-2RK
					L	21-U4-8LK	21-U4-2LK	22-U4-8LK	22-U4-2LK

下顎

部位	Torque	Angle	Rotation	M/D	R/L	.018	.022	.018	.022
前歯	-6	0	0	2.6mm	Uni	21-L1-8U	21-L1-2U	22-L1-8U	22-L1-2U
犬歯	-6	+3	0	3.4mm	R	21-L3-8R	21-L3-2R	22-L3-8R	22-L3-2R
					L	21-L3-8L	21-L3-2L	22-L3-8L	22-L3-2L
犬歯 フック付					R	21-L3-8RK	21-L3-2RK	22-L3-8RK	22-L3-2RK
					L	21-L3-8LK	21-L3-2LK	22-L3-8LK	22-L3-2LK
犬歯	0	+3	0	3.4mm	R	21-L3-8R0	21-L3-2R0	22-L3-8R0	22-L3-2R0
					L	21-L3-8L0	21-L3-2L0	22-L3-8L0	22-L3-2L0
犬歯 フック付					R	21-L3-8RK0	21-L3-2RK0	22-L3-8RK0	22-L3-2RK0
					L	21-L3-8LK0	21-L3-2LK0	22-L3-8LK0	22-L3-2LK0
第1小臼歯	-12	0	0	3.4mm	Uni	21-L4-8U	21-L4-2U	22-L4-8U	22-L4-2U
第1小臼歯 フック付					R	21-L4-8RK	21-L4-2RK	22-L4-8RK	22-L4-2RK
					L	21-L4-8LK	21-L4-2LK	22-L4-8LK	22-L4-2LK
第2小臼歯	-17	0	0	3.4mm	Uni	21-L5-8U	21-L5-2U	22-L5-8U	22-L5-2U
第2小臼歯 フック付					R	21-L5-8RK	21-L5-2RK	22-L5-8RK	22-L5-2RK
					L	21-L5-8LK	21-L5-2LK	22-L5-8LK	22-L5-2LK

症例セット(20個/セット)

内容	アンカーロックベース 定価6,500円		メッシュベース 定価7,900円	
	.018	.022	.018	.022
上下顎5~5 フックなし	21K-NHK-18	21K-NHK-22	22K-NHK-18	22K-NHK-22
上下顎5~5 3のみフック付	21K-3-18	21K-3-22	22K-3-18	22K-3-22
上下顎5~5 4, 5のみフック付	21K-45-18	21K-45-22	22K-45-18	22K-45-22
上下顎5~5 3, 4, 5フック付	21K-345-18	21K-345-22	22K-345-18	22K-345-22

* 全てのトルク、オフセット等の角度の単位は「°」です。




医療機器認証番号: 222AKBZX00123000

MAGNUMTM

BRACKET SERIES



by Orthodontic Design and Production, Inc.



Roth Prescription

上顎					10個/パック		10個/パック		10個/パック	
					スタンダードベース		オフセットベース		エクステンデッドベース	
部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
中切歯	+12	+5	3.9mm	R	30-U1-8R	30-U1-2R	---	---	---	---
				L	30-U1-8L	30-U1-2L	---	---	---	---
側切歯	+8	+9	3.0mm	R	30-U2-8R	30-U2-2R	---	---	---	---
				L	30-U2-8L	30-U2-2L	---	---	---	---
側切歯				R	30-U2-8RK	30-U2-2RK	---	---	---	---
フック付				L	30-U2-8LK	30-U2-2LK	---	---	---	---
犬歯	0	+11	3.5mm	R	30-U3-8R	30-U3-2R	---	---	---	---
				L	30-U3-8L	30-U3-2L	---	---	---	---
犬歯				R	30-U3-8RK	30-U3-2RK	---	---	---	---
フック付				L	30-U3-8LK	30-U3-2LK	---	---	---	---
小臼歯	-7	0	3.3mm	Uni	30-U4-8U	30-U4-2U	30-U4-8UO	30-U4-2UO	30-U4-8UE	30-U4-2UE
小臼歯				R	30-U4-8RK	30-U4-2RK	30-U4-8RKO	30-U4-2RKO	30-U4-8RKE	30-U4-2RKE
フック付				L	30-U4-8LK	30-U4-2LK	30-U4-8LKO	30-U4-2LKO	30-U4-8LKE	30-U4-2LKE

下顎

部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
前歯	0	0	2.5mm	Uni	30-L1-8U	30-L1-2U	---	---	---	---
犬歯	-11	+5	3.5mm	R	30-L3-8R	30-L3-2R	---	---	---	---
				L	30-L3-8L	30-L3-2L	---	---	---	---
犬歯				R	30-L3-8RK	30-L3-2RK	---	---	---	---
フック付				L	30-L3-8LK	30-L3-2LK	---	---	---	---
第1小臼歯	-17	0	3.3mm	Uni	30-L4-8U	30-L4-2U	30-L4-8UO	30-L4-2UO	30-L4-8UE	30-L4-2UE
第1小臼歯				R	30-L4-8RK	30-L4-2RK	30-L4-8RKO	30-L4-2RKO	30-L4-8RKE	30-L4-2RKE
フック付				L	30-L4-8LK	30-L4-2LK	30-L4-8LKO	30-L4-2LKO	30-L4-8LKE	30-L4-2LKE
第2小臼歯	-22	0	3.3mm	Uni	30-L5-8U	30-L5-2U	30-L5-8UO	30-L5-2UO	30-L5-8UE	30-L5-2UE
第2小臼歯				R	30-L5-8RK	30-L5-2RK	30-L5-8RKO	30-L5-2RKO	30-L5-8RKE	30-L5-2RKE
フック付				L	30-L5-8LK	30-L5-2LK	30-L5-8LKO	30-L5-2LKO	30-L5-8LKE	30-L5-2LKE

症例セット(20個/セット)

内容	スタンダードベース		オフセットベース		エクステンデッドベース	
	.018	.022	.018	.022	.018	.022
上下顎5~5 フックなし	30KR-NHK-18	30KR-NHK-22	30KO-NHK-18	30KO-NHK-22	30KE-NHK-18	30KE-NHK-22
上下顎5~5 3のみフック付	30KR-3-18	30KR-3-22	30KO-3-18	30KO-3-22	30KE-3-18	30KE-3-22
上下顎5~5 4, 5のみフック付	30KR-45-18	30KR-45-22	30KO-45-18	30KO-45-22	30KE-45-18	30KE-45-22
上下顎5~5 3, 4, 5フック付	30KR-345-18	30KR-345-22	30KO-345-18	30KO-345-22	30KE-345-18	30KE-345-22




* 全てのトルク、オフセット等の角度の単位は「°」です。

MAGNUMTM

BRACKET SERIES



by Orthodontic Design and Production, Inc.



Roth V Slot Prescription

					10個/パック		10個/パック		10個/パック	
上顎					スタンダードベース		オフセットベース		エクステンデッドベース	
部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
中切歯	+12	+5	3.9mm	R	31-U1-8R	31-U1-2R	---	---	---	---
				L	31-U1-8L	31-U1-2L	---	---	---	---
側切歯	+8	+9	3.0mm	R	31-U2-8R	31-U2-2R	---	---	---	---
				L	31-U2-8L	31-U2-2L	---	---	---	---
側切歯	0	+11	3.5mm	R	31-U2-8RK	31-U2-2RK	---	---	---	---
フック付				L	31-U2-8LK	31-U2-2LK	---	---	---	---
犬歯				R	31-U3-8R	31-U3-2R	---	---	---	---
				L	31-U3-8L	31-U3-2L	---	---	---	---
犬歯	0	+11	3.5mm	R	31-U3-8RK	31-U3-2RK	---	---	---	---
フック付				L	31-U3-8LK	31-U3-2LK	---	---	---	---
小臼歯	-7	0	3.3mm	Uni	31-U4-8U	31-U4-2U	31-U4-8UO	31-U4-2UO	31-U4-8UE	31-U4-2UE
小臼歯	-	0	3.3mm	R	31-U4-8RK	31-U4-2RK	31-U4-8RKO	31-U4-2RKO	31-U4-8RKE	31-U4-2RKE
フック付				L	31-U4-8LK	31-U4-2LK	31-U4-8LKO	31-U4-2LKO	31-U4-8LKE	31-U4-2LKE

下顎

部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
前歯	0	0	2.5mm	Uni	31-L1-8U	31-L1-2U	---	---	---	---
犬歯	-11	+5	3.5mm	R	31-L3-8R	31-L3-2R	---	---	---	---
				L	31-L3-8L	31-L3-2L	---	---	---	---
犬歯	-	5	3.5mm	R	31-L3-8RK	31-L3-2RK	---	---	---	---
フック付				L	31-L3-8LK	31-L3-2LK	---	---	---	---
第1小臼歯	-17	0	3.3mm	Uni	31-L4-8U	31-L4-2U	31-L4-8UO	31-L4-2UO	31-L4-8UE	31-L4-2UE
第1小臼歯	-	0	3.3mm	R	31-L4-8RK	31-L4-2RK	31-L4-8RKO	31-L4-2RKO	31-L4-8RKE	31-L4-2RKE
フック付				L	31-L4-8LK	31-L4-2LK	31-L4-8LKO	31-L4-2LKO	31-L4-8LKE	31-L4-2LKE
第2小臼歯	-22	0	3.3mm	Uni	31-L5-8U	31-L5-2U	31-L5-8UO	31-L5-2UO	31-L5-8UE	31-L5-2UE
第2小臼歯	-	0	3.3mm	R	31-L5-8RK	31-L5-2RK	31-L5-8RKO	31-L5-2RKO	31-L5-8RKE	31-L5-2RKE
フック付				L	31-L5-8LK	31-L5-2LK	31-L5-8LKO	31-L5-2LKO	31-L5-8LKE	31-L5-2LKE

症例セット(20個/セット)

内容	スタンダードベース		オフセットベース		エクステンデッドベース	
	.018	.022	.018	.022	.018	.022
上下顎5~5 フックなし	31KR-NHK-18	31KR-NHK-22	31KO-NHK-18	31KO-NHK-22	31KE-NHK-18	31KE-NHK-22
上下顎5~5 3のみフック付	31KR-3-18	31KR-3-22	31KO-3-18	31KO-3-22	31KE-3-18	31KE-3-22
上下顎5~5 4, 5のみフック付	31KR-45-18	31KR-45-22	31KO-45-18	31KO-45-22	31KE-45-18	31KE-45-22
上下顎5~5 3, 4, 5フック付	31KR-345-18	31KR-345-22	31KO-345-18	31KO-345-22	31KE-345-18	31KE-345-22


* 全てのトルク、オフセット等の角度の単位は"°"です。

MAGNUMTM

BRACKET SERIES



by Orthodontic Design and Production, Inc.



Hilgers Prescription

上顎					10個/パック スタンダードベース		10個/パック オフセットベース		10個/パック エクステンデッドベース	
部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
中切歯	+22	+5	3.9mm	R	33-U1-8R	33-U1-2R	---	---	---	---
				L	33-U1-8L	33-U1-2L	---	---	---	---
側切歯	+14	+9	3.0mm	R	33-U2-8R	33-U2-2R	---	---	---	---
				L	33-U2-8L	33-U2-2L	---	---	---	---
犬歯	7	+10	3.5mm	R	33-U3-8R	33-U3-2R	---	---	---	---
				L	33-U3-8L	33-U3-2L	---	---	---	---
犬歯 フック付				R	33-U3-8RK	33-U3-2RK	---	---	---	---
				L	33-U3-8LK	33-U3-2LK	---	---	---	---
小臼歯	-7	0	3.3mm	Uni	30-U4-8U	30-U4-2U	30-U4-8UO	30-U4-2UO	30-U4-8UE	30-U4-2UE
小臼歯				R	30-U4-8RK	30-U4-2RK	30-U4-8RKO	30-U4-2RKO	30-U4-8RKE	30-U4-2RKE
小臼歯 フック付				L	30-U4-8LK	30-U4-2LK	30-U4-8LKO	30-U4-2LKO	30-U4-8LKE	30-U4-2LKE

下顎

部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
前歯	0	0	2.5mm	Uni	30-L1-8U	30-L1-2U	---	---	---	---
犬歯	-7	+5	3.5mm	R	33-L3-8R	33-L3-2R	---	---	---	---
				L	33-L3-8L	33-L3-2L	---	---	---	---
犬歯 フック付				R	33-L3-8RK	33-L3-2RK	---	---	---	---
				L	33-L3-8LK	33-L3-2LK	---	---	---	---
第1小臼歯	-17	0	3.3mm	Uni	30-L4-8U	30-L4-2U	30-L4-8UO	30-L4-2UO	30-L4-8UE	30-L4-2UE
第1小臼歯				R	30-L4-8RK	30-L4-2RK	30-L4-8RKO	30-L4-2RKO	30-L4-8RKE	30-L4-2RKE
第1小臼歯 フック付				L	30-L4-8LK	30-L4-2LK	30-L4-8LKO	30-L4-2LKO	30-L4-8LKE	30-L4-2LKE
第2小臼歯	-22	0	3.3mm	Uni	30-L5-8U	30-L5-2U	30-L5-8UO	30-L5-2UO	30-L5-8UE	30-L5-2UE
第2小臼歯				R	30-L5-8RK	30-L5-2RK	30-L5-8RKO	30-L5-2RKO	30-L5-8RKE	30-L5-2RKE
第2小臼歯 フック付				L	30-L5-8LK	30-L5-2LK	30-L5-8LKO	30-L5-2LKO	30-L5-8LKE	30-L5-2LKE

症例セット(20個/セット)

内容	スタンダードベース		オフセットベース		エクステンデッドベース	
	.018	.022	.018	.022	.018	.022
上下顎5~5 フックなし	33KR-NHK-18	33KR-NHK-22	33KO-NHK-18	33KO-NHK-22	33KE-NHK-18	33KE-NHK-22
上下顎5~5 3のみフック付	33KR-3-18	33KR-3-22	33KO-3-18	33KO-3-22	33KE-3-18	33KE-3-22
上下顎5~5 4, 5のみフック付	33KR-45-18	33KR-45-22	33KO-45-18	33KO-45-22	33KE-45-18	33KE-45-22
上下顎5~5 3, 4, 5フック付	33KR-345-18	33KR-345-22	33KO-345-18	33KO-345-22	33KE-345-18	33KE-345-22

* 全てのトルク、オフセット等の角度の単位は「°」です。




医療機器認証番号: 222AKBZX00123000

MAGNUMTM

BRACKET SERIES



by Orthodontic Design and Production, Inc.



Ricketts Prescription

					10個/パック		10個/パック		10個/パック	
上顎					スタンダードベース		オフセットベース		エクステンデッドベース	
部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
中切歯	+22	+5	3.9mm	R	33-U1-8R	33-U1-2R	---	---	---	---
				L	33-U1-8L	33-U1-2L	---	---	---	---
側切歯	+14	+9	3.0mm	R	33-U2-8R	33-U2-2R	---	---	---	---
				L	33-U2-8L	33-U2-2L	---	---	---	---
犬歯	+7	+5	3.5mm	R	34-U3-8R	34-U3-2R	---	---	---	---
				L	34-U3-8L	34-U3-2L	---	---	---	---
犬歯				R	34-U3-8RK	34-U3-2RK	---	---	---	---
フック付				L	34-U3-8LK	34-U3-2LK	---	---	---	---
小臼歯	0	0	3.3mm	Uni	34-U4-8U	34-U4-2U	34-U4-8UO	34-U4-2UO	34-U4-8UE	34-U4-2UE
小臼歯				R	34-U4-8RK	34-U4-2RK	34-U4-8RKO	34-U4-2RKO	34-U4-8RKE	34-U4-2RKE
フック付				L	34-U4-8LK	34-U4-2LK	34-U4-8LKO	34-U4-2LKO	34-U4-8LKE	34-U4-2LKE

下顎

部位	Torque	Angle	M/D	R/L	.018	.022	.018	.022	.018	.022
前歯	0	0	2.5mm	Uni	30-L1-8U	30-L1-2U	---	---	---	---
犬歯	-7	+5	3.5mm	R	33-L3-8R	33-L3-2R	---	---	---	---
				L	33-L3-8L	33-L3-2L	---	---	---	---
犬歯				R	33-L3-8RK	33-L3-2RK	---	---	---	---
フック付				L	33-L3-8LK	33-L3-2LK	---	---	---	---
第1・第2小臼歯	0	0	3.3mm	Uni	34-L4-8U	34-L4-2U	34-L4-8UO	34-L4-2UO	34-L4-8UE	34-L4-2UE
第1・第2小臼歯				R	34-L4-8RK	34-L4-2RK	34-L4-8RKO	34-L4-2RKO	34-L4-8RKE	34-L4-2RKE
フック付				L	34-L4-8LK	34-L4-2LK	34-L4-8LKO	34-L4-2LKO	34-L4-8LKE	34-L4-2LKE

フルセット(20個入)

内容	スタンダードベース		オフセットベース		エクステンデッドベース	
	.018	.022	.018	.022	.018	.022
上下顎5~5 フックなし	34KR-NHK-18	34KR-NHK-22	34KO-NHK-18	34KO-NHK-22	34KE-NHK-18	34KE-NHK-22
上下顎5~5 3のみフック付	34KR-3-18	34KR-3-22	34KO-3-18	34KO-3-22	34KE-3-18	34KE-3-22
上下顎5~5 4, 5のみフック付	34KR-45-18	34KR-45-22	34KO-45-18	34KO-45-22	34KE-45-18	34KE-45-22
上下顎5~5 3, 4, 5フック付	34KR-345-18	34KR-345-22	34KO-345-18	34KO-345-22	34KE-345-18	34KE-345-22


MAGNUMTM

BRACKET SERIES



by Orthodontic Design and Production, Inc.

ウェルドタイプ 小臼歯ブラケット



上顎

10個/パック

部位	Torque	Angle	M/D	R/L	.018	.022
小臼歯	0	0	3.3mm	Uni	34-WU4-8U	34-WU4-2U
	-7	0	3.3mm	Uni	30-WU4-8U	30-WU4-2U
小臼歯	0	0	3.3mm	R	34-WU4-8RK	34-WU4-2RK
フック付				L	34-WU4-8LK	34-WU4-2LK
	-7	0	3.3mm	R	30-WU4-8RK	30-WU4-2RK
				L	30-WU4-8LK	30-WU4-2LK

下顎

部位	Torque	Angle	M/D	R/L	.018	.022
第1小臼歯	0	0	3.3mm	Uni	34-WL4-8U	34-WL4-2U
	-17	0	3.3mm	Uni	30-WL4-8U	30-WL4-2U
第1小臼歯	0	0	3.3mm	R	34-WL4-8RK	34-WL4-2RK
フック付				L	34-WL4-8LK	34-WL4-2LK
	-17	0	3.3mm	R	30-WL4-8RK	30-WL4-2RK
				L	30-WL4-8LK	30-WL4-2LK
第2小臼歯	0	0	3.3mm	Uni	34-WL4-8U	34-WL4-2U
	-22	0	3.3mm	Uni	30-WL5-8U	30-WL5-2U
第2小臼歯	0	0	3.3mm	R	34-WL4-8RK	34-WL4-2RK
フック付				L	34-WL4-8LK	34-WL4-2LK
	-22	0	3.3mm	R	30-WL5-8RK	30-WL5-2RK
				L	30-WL5-8LK	30-WL5-2LK

バーチカルスロットタイプ



10個/パック

部位	Torque	Angle	M/D	R/L	.018	.022
上顎小臼歯	-7	0	3.3mm	Uni	31-WU4-8U	31-WU4-2U
下顎小臼歯	-17	0	3.3mm	Uni	31-WL4-8U	31-WL4-2U

* 全てのトルク、オフセット等の角度の単位は「°」です。

医療機器認証番号: 222AKBZX00123000